

Ideal For Home & Travel

Lightweight & Portable!

Includes dual purpose cleaning and travel case for at home or on the go sanitizing!

Easy To Use

Long lasting rechargeable lithium-ion battery with USB charging adapter.

Maintenance Free

No additional pre-wash solutions, filters, or replacement bags required!

Purify O3 Makes Sanitizing Easy!



Attach mask & tubing to Purify O3. (Surface clean with damp cloth to remove residue.)



Place mask/tubing assembly and Purify O3 in case. Press and hold button to start, zip case shut.



After 35 minute cycle, unit will auto shut off. Let rest for additional 90 minutes to allow ozone to dissipate.



833-825-0672 www.shoprotech.com

2 Year Warranty - Register Online At: respondo2.com/warranty

Purify Q3

Portable Ozone CPAP/BiPAP Sanitizer





Sanitize with Purify Q

The simple, effective way to clean your CPAP/BiPAP mask & equipment

www.purifyo3.com

REST EASY WITH



Disinfect With Just The Touch Of A Button

Simple, Easy To Use Controls: One button push to start with automatic shut-off.

Quiet, Worry-Free Cleaning: Whisper quiet design, audible alert when cycle is complete.

Universal Design: Fits all 22mm CPAP/BiPAP equipment and masks.

Safe & Effective: Provides hospital grade sanitation in just 35 minutes. FDA Registered.

All-In-One Solution

Patient friendly sanitizing kit includes Purify 03, storage/sanitizing case, charging cable, and tubing adapter.





Why Sanitize With Ozone Technology?

Natural: Activated oxygen, known as ozone (O3), is a naturally formed gas and is widely used as a disinfecting agent in a number of industries world-wide.

Effective: Purify O3 provides hospital grade sanitation and quickly, permanently kills 99.99% of viruses and bacteria.

Affordable: No pre-wash, replacement bags, or extra filters required. Purify O3 includes all you need to sanitize your device in one affordable kit! HSA/FSA Eligible.